

Super Cut Scissors



Some of the benefits of Super-Cut Scissors:

- Cut through soft and fibrous tissue without effort.
- Superb control permits precision in dissection.
- Less hand fatigue for the surgeon and less trauma for the patient.
- Retains sharpness longer than ordinary scissors.
- Shafts are fully gold plated for easy identification.
- Fully autoclavable.
- Not recommended for fiberglass or plaster splints.



PROFESSIONAL Supercut scissors with their uniquely ground edges have with ultimate cutting efficiency. The special micro serrations in conjunction with the uniquely ground edges ensure precise, fine and minimally traumatising cuts. Professional supercut scissors have fully gold plated shafts for ease of identification.

The range of supercut scissors from PROFESSIONAL comprises all models used routinely by surgeons: scissors for delicate and extremely delicate dissecting, e.g. Iris and Micro, models according to Metzenbaum in all possible lengths and designs, e.g. in a slim version for especially fine cutting to name but a few.

The supercut scissors can also be supplied with "TC" welded cutting edges of solid tungsten carbide alloy. The junction between the steel of the scissors and the tungsten carbide is seamless and homogenous. The welded cutting edges of solid tungsten carbide alloy guarantee excellent cutting quality and long life-span.

Identification:

PROFESSIONAL TC scissors:
gold-coloured rings



PROFESSIONAL supercut scissors:
one shaft completely gold plated




PROFESSIONAL supercut plus line scissors:
one shaft completely gold plated and
one ring gold plated



Supercut Scissors

Iris Scissors



Straight			
	10.5 cm/4 $\frac{1}{8}$ "	11.5 cm/4 $\frac{1}{2}$ "	12.5 cm/5"
TC		PC-05.61	
SUPERCUT	PC-01.110	PC-01.111	
SUPERCUT PLUS	PC-01.210	PC-05.161	

Curved			
TC		PC-06.61	PC-02.112
SUPERCUT	PC-02.110	PC-02.111	
SUPERCUT PLUS	PC-02.210	PC-06.161	


IRIS SCISSORS

Popular sharp/sharp scissors offer outstanding distal cutting capability of tissue and suture.

Note: SuperCut scissors are not recommended for use on suture material.




MICRO IRIS
PC-04.130
9 cm/3 $\frac{1}{2}$ "
SUPERCUT



Straight	
	11 cm/4 $\frac{3}{8}$ "
SUPERCUT	PC-64.111
BLACK LINE	PC-64.210

Curved	
SUPERCUT	PC-64.112
BLACK LINE	PC-64.212



Straight	
	15 cm/6"
SUPERCUT	PC-65.114
BLACK LINE	PC-65.214

Curved	
SUPERCUT	PC-66.114
BLACK LINE	PC-66.214

STEVENS TENOTOMY SCISSORS

The delicate blade design tapers for fine tissue penetration and dissection. Tenotomy scissors are excellent for general, facial, ophthalmic, dermatological and hand surgery procedures.



STEVENS


Straight	11.5 cm/4½"
SUPERCUT	PC-10.112

Curved	
SUPERCUT	PC-10.113



STEVENS


Straight	13 cm/5½"
SUPERCUT	PC-10.115



KNAPP

Straight	12.5 cm/5"
TC	PC-13.62
SUPERCUT	PC-11.112
SUPERCUT PLUS	PC-13.162

Curved	
TC	PC-13.63
SUPERCUT	PC-12.112
SUPERCUT PLUS	PC-13.163



IRIS
RIBBON HANDLE
SUPERCUT

9.5 cm/ 3¾"	10.5 cm/ 4¼"	9.5 cm/ 3¾"	10.5 cm/ 4¼"
PC-07.109	PC-07.110	PC-08.109	PC-08.110
Straight	Straight	Curved	Curved

Supercut Scissors

Operating Scissors

Blunt / Blunt




Straight		
	13 cm/4 $\frac{1}{8}$ "	14.5 cm/5 $\frac{3}{4}$ "
TC		PC-14.64
SUPERCUT	PC-14.113	PC-14.114
Curved		
TC		PC-18.64
SUPERCUT	PC-18.113	PC-18.114

Sharp / Sharp



Straight		
	13 cm/4 $\frac{1}{8}$ "	14.5 cm/5 $\frac{3}{4}$ "
TC		PC-15.64
SUPERCUT	PC-15.113	PC-15.114
Curved		
TC		PC-19.64
SUPERCUT	PC-19.113	PC-19.114

Sharp / Blunt




Straight		
	13 cm/4 $\frac{1}{8}$ "	14.5 cm/5 $\frac{3}{4}$ "
TC		PC-16.64
SUPERCUT	PC-16.113	PC-16.114
Curved		
TC		PC-17.64
SUPERCUT	PC-17.113	PC-17.114

General utility scissors with sharp-blunt and blunt-blunt tips are ideal for cutting suture, gauze, tubing, etc. The sharp-sharp patterns are also used for general utility applications but preferred for cutting tendons and cartilage.

Note: SuperCut scissors are not recommended for use on suture, gauze or tubing.

Dissecting scissors




Sturdy, durable and reliable scissors. The robust and heavier design is used to precisely cut and dissect tough, thick and muscular tissue.

MAYO

Straight	14 cm/5½"	15.5 cm/6¼"	17 cm/6¾"	23 cm/9"
TC	PC-24.64			
SUPERCUT	PC-24.114	PC-24.115	PC-24.117	PC-24.123
SUPERCUT PLUS	PC-24.164	PC-24.165	PC-24.167	PC-24.173

Curved				
TC	PC-25.64			
SUPERCUT	PC-25.114	PC-25.115	PC-25.117	PC-25.123
SUPERCUT PLUS	PC-25.164	PC-25.165	PC-25.167	PC-25.173

reinforced rounded blades provide additional strength and durability when used to cut or dissect the very toughest tissue.



MAYO-STILLE

Straight	15 cm/6"	17 cm/6¾"
SUPERCUT	PC-28.115	PC-28.117
BLACK LINE	PC-28.215	PC-28.217
SUPERCUT PLUS	PC-29.165	PC-29.167


Curved		
SUPERCUT	PC-30.115	PC-30.117
BLACK LINE	PC-30.215	PC-30.217
SUPERCUT PLUS	PC-31.165	PC-31.167

Supercut Scissors

Abdominal Scissors




DOYEN




Straight		18 cm/7 ¹ / ₈ "
SUPERCUT		PC-22.118
BLACK LINE		PC-22.218
Curved		
SUPERCUT		PC-23.118
BLACK LINE		PC-23.218



SIMS




Straight		20 cm/8"	23 cm/9"
SUPERCUT		PC-26.120	
BLACK LINE			
SUPERCUT PLUS			
Curved			
SUPERCUT		PC-27.120	PC-27.123
BLACK LINE		PC-27.220	PC-27.223
SUPERCUT PLUS		PC-27.150	



SIMS-SIEBOLD
PC-27.125
25 cm/10"
SUPERCUT




SIEBOLD
PC-27.124
24.5 cm/9¾"
SUPERCUT




RESANO
25 cm/10"
SUPERCUT PC-27.135
BLACK LINE PC-27.235

Supercut Scissors

Dissecting Scissors




MAYO




18 cm/7 1/8"
PC-29.118


SUPERCUT



23 cm/9"
PC-29.123




MAYO



18 cm/7 1/8"
PC-29.138

SUPERCUT



23 cm/9"
PC-29.143

Dissecting scissors



MAYO-LEXER




Straight	16 cm/6¼"
SUPERCUT	PC-32.116
BLACK LINE	PC-32.216
SUPERCUT PLUS	PC-34.196

Curved	
SUPERCUT	PC-32.117
BLACK LINE	PC-32.217
SUPERCUT PLUS	PC-35.196




LEXER



	16.5 cm/ 6½"
SUPERCUT	PC-32.187
BLACK LINE	PC-32.287
SUPERCUT PLUS	PC-32.197



MINI-METZENBAUM



Straight	9 cm/3½"
SUPERCUT	PC-36.109
BLACK LINE	PC-36.209

Curved	
SUPERCUT	PC-37.109
BLACK LINE	PC-37.209



BABY-METZENBAUM



Straight	11.5 cm/4½"
SUPERCUT	PC-36.111
BLACK LINE	PC-36.211
SUPERCUT PLUS	PC-36.161

Curved	
SUPERCUT	PC-37.111
BLACK LINE	PC-37.211
SUPERCUT PLUS	PC-37.161


Supercut Scissors

Dissecting Scissors



METZENBAUM

This very popular and versatile pattern features gently rounded tips for dissection. a light and perfectly balanced design precisely cuts either tough or delicate tissue.




	Straight	14 cm/5½"	15 cm/6"	18 cm/7"	20 cm/8"
SUPERCUT	PC-36.114			PC-38.118	PC-38.120
BLACK LINE	PC-36.214			PC-38.218	PC-38.220
SUPERCUT PLUS	PC-36.164			PC-38.168	PC-38.170

	Curved				
SUPERCUT	PC-37.114	PC-37.115	PC-39.118	PC-39.120	
BLACK LINE	PC-37.214	PC-37.215	PC-39.218	PC-39.220	
SUPERCUT PLUS	PC-37.164		PC-39.168	PC-39.170	



NELSON-METZENBAUM




	Straight	23 cm/8"	26 cm/10¼"	30 cm/11¾"
SUPERCUT	PC-40.123		PC-40.126	PC-40.130
BLACK LINE	PC-40.223		PC-40.226	PC-40.230
SUPERCUT PLUS	PC-40.173			

	Curved			
SUPERCUT	PC-41.123	PC-41.126	PC-41.130	
BLACK LINE	PC-41.223	PC-41.226	PC-41.230	
SUPERCUT PLUS	PC-41.173	PC-41.176	PC-41.180	



METZENBAUM-SLIM



Straight	14 cm/5½"	18 cm/7"	20 cm/8"	23 cm/9"
SUPERCUT	PC-42.114	PC-44.118	PC-44.120	PC-44.123
BLACK LINE	PC-42.214	PC-44.218	PC-44.220	PC-44.223
SUPERCUT PLUS	PC-42.164	PC-44.168	PC-44.170	PC-44.173

Curved				
SUPERCUT	PC-43.114	PC-45.118	PC-45.120	PC-45.123
BLACK LINE	PC-43.214	PC-45.218	PC-45.220	PC-45.223
SUPERCUT PLUS	PC-43.164	PC-45.168	PC-45.170	PC-45.173




METZENBAUM-SLIM




	26 cm/10¼"	28.5 cm/11¼"
SUPERCUT	PC-45.176	PC-45.178

Supercut Scissors


Dissecting Scissors




METZENBAUM




	18 cm/7 $\frac{1}{8}$ "	20 cm/8"	23 cm/9"
SUPERCUT	PC-46.168	PC-46.170	PC-46.173
BLACK LINE	PC-46.268	PC-46.270	PC-46.273
SUPERCUT PLUS	PC-46.178	PC-46.180	PC-46.183



METZENBAUM-FREEMAN




	18 cm/7 $\frac{1}{8}$ "	20 cm/8"	23 cm/9"
SUPERCUT	PC-47.168	PC-47.170	PC-47.123
BLACK LINE	PC-47.268	PC-47.270	PC-47.273




KELLY

Straight	16 cm/6¼"
SUPERCUT	PC-46.116
BLACK LINE	PC-46.216
SUPERCUT PLUS	PC-46.166
Curved	
SUPERCUT	PC-47.116
BLACK LINE	PC-47.216
SUPERCUT PLUS	PC-47.166




BRÜSER-KELLY

Straight	16 cm/6¼"
SUPERCUT	PC-48.116
Curved	
SUPERCUT	PC-49.116



TOENNIS-ADSON

	18cm/7½"
SUPERCUT	PC-50.118
BLACK LINE	PC-50.218
SUPERCUT PLUS	PC-50.168



	15cm/6"
SUPERCUT	PC-50.125
BLACK LINE	PC-50.225



	15cm/6"
SUPERCUT	PC-51.115
BLACK LINE	PC-51.215

Supercut Scissors

Dissecting Scissors




	16.5cm/6½"
SUPERCUT	PC-57.116
BLACK LINE	PC-57.216



SIEBERT

Straight	11.5 cm/4½"
SUPERCUT	PC-58.131

Curved	
SUPERCUT	PC-59.132




RAGNELL (KILNER)

	13 cm/5½"	15 cm/6"	18 cm/7½"
SUPERCUT	PC-59.113	PC-59.115	PC-59.118
BLACK LINE	PC-59.213	PC-59.215	PC-59.218
SUPERCUT PLUS	PC-59.163	PC-59.165	PC-59.168



GREENBERG PAR


	11.5 cm/4½"
SUPERCUT	PC-60.122
BLACK LINE	PC-60.222
SUPERCUT PLUS	PC-61.122



PECK-JOSEPH

Straight	14 cm/5½"
SUPERCUT	PC-65.134
BLACK LINE	PC-65.234


Curved	
SUPERCUT	PC-66.134
BLACK LINE	PC-66.234



PECK-JOSEPH


Straight	14 cm/5½"
SUPERCUT	PC-65.144
BLACK LINE	PC-65.244

Curved	
SUPERCUT	PC-66.144
BLACK LINE	PC-66.244




Straight	14 cm/5½"	18 cm/7½"
SUPERCUT	PC-67.114	
SUPERCUT PLUS	PC-67.164	

Curved		
SUPERCUT	PC-68.114	PC-68.118
SUPERCUT PLUS	PC-68.164	



JAMESON

	15 cm/6"
SUPERCUT	PC-69.117
BLACK LINE	PC-69.217



JAMESON


	18 cm/7½"
SUPERCUT	PC-70.118
BLACK LINE	PC-70.218

Supercut Scissors

Dissecting Scissors




REYNOLDS (JAMESON)



	12 cm/4¾"	16 cm/6¼"	18 cm/7½"
SUPERCUT	PC-73.112	PC-73.116	PC-73.118
BLACK LINE	PC-73.212	PC-73.216	PC-73.218




FANOUS




Straight	16.5 cm/6½"	20.5 cm/8"
SUPERCUT	PC-74.116	PC-74.126

Curved		
SUPERCUT	PC-74.120	PC-74.130

Dissecting scissors




GORNEY
FACELIFT




	17 cm/6¾"	19 cm/7½"	23 cm/9"
SUPERCUT	PC-77.127	PC-77.129	PC-77.133
BLACK LINE	PC-77.227	PC-77.229	PC-77.233

The Gorney & Gorney-freeman scissors are a heavier pattern when compared to the Kaye and Kaye-freeman scissors and are ideal for face lift and general surgery procedures. The saber-back blades with semi-sharp outer edges combined with the flattened tip allow for effortlessly passing through and undermining tissue on the face, neck, etc. The flattened and delicate blunt tips facilitate precision dissection and eliminate potential for "buttonholing".




GORNEY-FREEMAN
FACELIFT



	17 cm/6¾"	19 cm/7½"	23 cm/9"
SUPERCUT	PC-80.117	PC-80.119	PC-80.123
BLACK LINE	PC-80.217	PC-80.219	PC-80.223


Supercut Scissors

Dissecting Scissors



	23 cm/9"
SUPERCUT	PC-80.152
BLACK LINE	PC-80.252

GORNEY-FREEMAN



	11.5 cm/4½"
SUPERCUT	PC-81.132
BLACK LINE	PC-81.232


KAYE



Straight	13 cm/5½"
SUPERCUT	PC-90.113
BLACK LINE	PC-90.213
SUPERCUT PLUS	PC-90.163

FOX

Curved	
SUPERCUT	PC-91.113
BLACK LINE	PC-91.213
SUPERCUT PLUS	PC-91.163



	15 cm/6"	17 cm/6¾"	19 cm/7½"	23 cm/9"
SUPERCUT	PC-81.115	PC-81.117	PC-81.119	PC-81.123
BLACK LINE	PC-81.215	PC-81.217	PC-81.219	PC-81.223

NUMERICAL INDEX

PC-01.110	2	PC-24.217	5	PC-37.109	9	PC-45.168	11	PC-64.212	2	PC-81.223	18
PC-01.111	2	PC-24.223	5	PC-37.111	9	PC-45.170	11	PC-65.114	2	PC-81.232	18
PC-01.210	2	PC-25.114	5	PC-37.114	10	PC-45.173	11	PC-65.134	15	PC-90.113	18
PC-01.211	2	PC-25.115	5	PC-37.115	10	PC-45.176	11	PC-65.144	15	PC-90.163	18
PC-02.110	2	PC-25.117	5	PC-37.161	9	PC-45.178	11	PC-65.214	2	PC-90.213	18
PC-02.111	2	PC-25.123	5	PC-37.164	10	PC-45.218	11	PC-65.234	15	PC-91.113	18
PC-02.210	2	PC-25.164	5	PC-37.209	9	PC-45.220	11	PC-65.244	15	PC-91.163	18
PC-02.211	2	PC-25.165	5	PC-37.211	9	PC-45.223	11	PC-66.114	2	PC-91.213	18
PC-04.130	2	PC-25.167	5	PC-37.214	10	PC-46.116	13	PC-66.134	15		
PC-05.161	2	PC-25.173	5	PC-37.215	10	PC-46.166	13	PC-66.144	15		
PC-06.161	2	PC-25.214	5	PC-38.118	10	PC-46.168	12	PC-66.214	2		
PC-10.112	3	PC-25.217	5	PC-38.120	10	PC-46.170	12	PC-66.234	15		
PC-10.113	3	PC-25.223	5	PC-38.168	10	PC-46.173	12	PC-66.244	15		
PC-10.115	3	PC-26.120	6	PC-38.170	10	PC-46.178	12	PC-67.114	15		
PC-10.115	3	PC-27.120	6	PC-38.218	10	PC-46.180	12	PC-67.164	15		
PC-10.212	3	PC-27.123	6	PC-38.220	10	PC-46.183	12	PC-67.164	15		
PC-10.213	3	PC-27.124	7	PC-39.118	10	PC-46.216	13	PC-68.114	15		
PC-10.215	3	PC-27.125	7	PC-39.120	10	PC-46.268	12	PC-68.114	15		
PC-10.215	3	PC-27.135	7	PC-39.168	10	PC-46.270	12	PC-68.118	15		
PC-11.112	3	PC-27.150	6	PC-39.170	10	PC-46.273	12	PC-68.164	15		
PC-11.212	3	PC-27.220	6	PC-39.218	10	PC-47.116	13	PC-69.117	15		
PC-12.112	3	PC-27.223	6	PC-39.220	10	PC-47.123	12	PC-69.217	15		
PC-12.212	3	PC-27.235	7	PC-40.123	10	PC-47.166	13	PC-70.118	15		
PC-13.162	3	PC-28.115	5	PC-40.126	10	PC-47.168	12	PC-70.218	15		
PC-13.163	3	PC-28.117	5	PC-40.130	10	PC-47.170	12	PC-73.112	16		
PC-14.113	4	PC-28.215	5	PC-40.173	10	PC-47.216	13	PC-73.116	16		
PC-14.114	4	PC-28.217	5	PC-40.223	10	PC-47.268	12	PC-73.118	16		
PC-14.214	4	PC-29.118	8	PC-40.226	10	PC-47.270	12	PC-73.212	16		
PC-15.113	4	PC-29.123	8	PC-40.230	10	PC-47.273	12	PC-73.216	16		
PC-15.114	4	PC-29.138	8	PC-41.123	10	PC-48.116	13	PC-73.218	16		
PC-15.214	4	PC-29.143	8	PC-41.126	10	PC-49.116	13	PC-74.116	16		
PC-16.113	4	PC-29.165	5	PC-41.130	10	PC-50.118	13	PC-74.120	16		
PC-16.114	4	PC-29.167	5	PC-41.173	10	PC-50.125	13	PC-74.126	16		
PC-16.214	4	PC-30.115	5	PC-41.176	10	PC-50.168	13	PC-74.130	16		
PC-17.113	4	PC-30.117	5	PC-41.180	10	PC-50.218	13	PC-77.127	17		
PC-17.114	4	PC-30.215	5	PC-41.223	10	PC-50.225	13	PC-77.129	17		
PC-17.214	4	PC-30.217	5	PC-41.226	10	PC-51.115	13	PC-77.133	17		
PC-18.113	4	PC-31.165	5	PC-41.230	10	PC-51.215	13	PC-77.227	17		
PC-18.114	4	PC-31.167	5	PC-42.114	11	PC-57.116	14	PC-77.229	17		
PC-18.214	4	PC-32.116	9	PC-42.164	11	PC-57.216	14	PC-77.233	17		
PC-19.113	4	PC-32.117	9	PC-42.214	11	PC-58.131	14	PC-80.117	17		
PC-19.114	4	PC-32.187	9	PC-43.114	11	PC-59.113	14	PC-80.119	17		
PC-19.214	4	PC-32.197	9	PC-43.164	11	PC-59.115	14	PC-80.123	17		
PC-22.118	6	PC-32.216	9	PC-43.214	11	PC-59.118	14	PC-80.152	18		
PC-22.218	6	PC-32.217	9	PC-44.118	11	PC-59.163	14	PC-80.217	17		
PC-23.118	6	PC-32.287	9	PC-44.120	11	PC-59.165	14	PC-80.219	17		
PC-23.218	6	PC-34.196	9	PC-44.123	11	PC-59.168	14	PC-80.223	17		
PC-24.114	5	PC-35.196	9	PC-44.168	11	PC-59.213	14	PC-80.252	18		
PC-24.115	5	PC-36.109	9	PC-44.170	11	PC-59.215	14	PC-81.115	18		
PC-24.117	5	PC-36.111	9	PC-44.173	11	PC-59.218	14	PC-81.117	18		
PC-24.123	5	PC-36.114	10	PC-44.218	11	PC-60.122	14	PC-81.119	18		
PC-24.164	5	PC-36.161	9	PC-44.220	11	PC-60.222	14	PC-81.123	18		
PC-24.165	5	PC-36.164	10	PC-44.223	11	PC-61.122	14	PC-81.132	18		
PC-24.167	5	PC-36.209	9	PC45.118	11	PC-64.111	2	PC-81.215	18		
PC-24.173	5	PC-36.211	9	PC-45.120	11	PC-64.112	2	PC-81.217	18		
PC-24.214	5	PC-36.214	10	PC-45.123	11	PC-64.210	2	PC-81.219	18		



Address:
PROFESSIONAL HOSPITAL FURNISHERS
19-A (N), Industrial Estate,
P. O. Box: 2095, Sialkot, 51310-Pakistan.
Tel. +92 52 325 6440, 325 6441
Fax: +92 52 325 6442, 355 2939
Email : info@phfworld.com